OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands
T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

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OmniMedical
the WISH company

OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands

T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

Instructions for Use – Comfort Washable Underpads

For:
• Comfort washable underpad 70x95 cm
• Comfort Plus washable underpad with tuck-in flaps 80x90 cm

Application

The Comfort washable underpads are designed to protect bedding, mattresses, and chairs against moisture loss such as urine or other bodily fluids. Suitable for home care, long-term care, and private use.

Suitable for

• Light to heavy incontinence
• Heavy menstruation
• Maternity care and childbirth
• Use on beds, chairs, commodes, or in baby cribs
• Use at home or while travelling/on holiday

Before First Use

Always wash the underpad before first use. This improves absorbency and ensures optimal comfort.

Instructions for Use

1. Place the underpad with the soft side facing up.
2. Position the underpad on the desired surface.
3. For the Comfort Plus version, the tuck-in flaps can be secured under the mattress for extra stability.
4. Check the underpad regularly for saturation.
5. Replace or wash the underpad when necessary.

Washing Instructions

Washing

• Recommended washing temperature: 60°C – 90°C
• Maximum washable temperature: 90°C
• Suitable for intensive washing, both at home and professionally

OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands
T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

Drying

- Suitable for tumble drying
- Maximum drying temperature: 100°C

Ironing

- Iron only with a cold iron.

Do Not Use

- Do not use bleach or aggressive chemicals
- Do not dry clean

Product Features

Comfort Underpad 70x95 cm

- Two-layer design
- Absorbency: approx. 700 ml
- Soft, skin-friendly top layer
- Flexible and leak-proof backing material

Comfort Plus Underpad 80x90 cm

- Four-layer design
- Absorbency: approx. 1800 ml
- Reduces odour and fluid backflow
- Equipped with convenient tuck-in flaps
- Extra stable and comfortable

Safety Information

- Regularly inspect the underpad for wear or damage.
- Use only as intended according to the specified application.

Storage

- Store in a dry and clean place
- Avoid prolonged exposure to direct sunlight