OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands
T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

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Instruction for use Comfort Urinal

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1.2 IFU (Lay Person/Home Use)

Before use	-	Check the urinal for damage before using it.
	-	Do not use an unsuitable or damaged urinal!
How to use	-	Hold or place the urinal between your legs.
	-	Check that your urine cannot flow over the edge of the urinal.
After use	-	Empty the urinal into the toilet or the toilet flusher.
	-	Clean (or wash) the urinal after use using disinfectants in the Ctgb database (Appendix 2).
	-	Place the urinal back in the holder.
Cross-contamination	-	The urinal is multiple use for one patient only. Do not share the urinal with other persons to protect yourself against potential cross-infections.
Storage	-	Dry storage between 5-40°Celsius.
	-	Avoid direct sun exposure during storage.
Shelf life	-	Based on comparable urinals made of PP and PC, the urinal has a shelf life of approximately five (5) years.
Disposal	-	Follow applicable (local) regulations for disposal.

OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands
T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

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2.2 IFU (Professional User)

Indications	- The Comfort (urine bottle) enables bedridden people to pass urine while lying down.
Before use	- Checks the urinal for damage before using it. - Do not use an unsuitable or damaged urinal.
How to use	- Empty the urinal into the toilet or the toilet flusher, or use the internal process to empty the collected urine.
Single-patient	- The urinal is multiple use for one patient only. Do not use the urinal for other persons.
Cross-contamination	- A potential risk is cross-contamination or transmission of pathogenic microorganisms when multiple patients use.
Washer-disinfector	- Use a washer-disinfector to clean and disinfect the urinal before use. - The washer-disinfector must comply with NEN-EN-ISO 15883-3:2006 as thermal disinfection. In this standard are applicable requirements and criteria specified.
Thermal disinfection	- When using thermal disinfection of the urinal. When professional use: - read WIP: 161220 RL Urinelozing-stoelgang_wijziging dec16DEF - read WIP: 161220 Pospoelers Generiek_DEFINITIEF
Sterilization	- The urine bottles (PP/PC) can be autoclaved (heat sterilization) up to 121°C. - Read the autoclave manual before use and follow the instructions belonging to the autoclave. - When using an autoclave, the PP/PC surface may change due to thermal influences. - Sterilization temperature: 121°C (chamber pressure approx. 1.1 bar). (duration: 15 min). - The autoclave cleaning duration depends on the manufacturer's instructions and is also based on the recommended cleaning agents and disinfectants.
Storage	- Dry storage between 5-40°C - Avoid direct sun exposure.
Shelf life	- Based on comparable urinals made of PP and PC: the urinal has a shelf life of approximately five (5) years.
Disposal	- Products contaminated by chemical or (human) biological substances must be disposed of by following (local) applicable regulations for the relevant chemicals and/or infectious material.

OmniMedical BV
Weverstraat 28
4061 AP Ophemert
The Netherlands
T +31 (0)344 65 14 44
info@omnimedical.nl
www.omnimedical.nl

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Appendix 2

Search strategy overview approved disinfectants

- Go to www.ctgb.nl.
- Under the heading 'Authorisation', select 'Download permitted and expired substances'.
- Click on 'Download overview of authorized substances'. An Excel file will now open.
- This Excel file contains all permitted resources. Important information in this Excel file:
 In column B you will find the name of the resources.
 Column C contains links to the current instructions for use of the disinfectants. Copy a link into your internet browser to open the file.
 In column T you will find which PT code(s) applies to the product. The PT code indicates for which material the product is suitable. Agents suitable for hand disinfection have a PT01 code (''Biocides for human hygiene'). Agents suitable for materials and surfaces have a PT02 code 'Disinfectants for private use and for public health, as well as other disinfectants'.